

BPL120/80_{B1}



Happier Living Everyday

Count it Easy



Pressure Check Made Easy



BPL 120/80_{B1}



Clinically accurate blood pressure readings made possible with easy operations

Features



Automatic Pressurization

High-tech measurement technology to monitor blood pressure from the arm – all with the push of a button



Large LCD Display

Extra large LCD display of 84.1 X 55.1 mm provides excellent display clarity and enables easy reading



Irregular Heartbeat Detection

Automatically alerts the user when heart rhythm is either 25% faster or slower than the average rhythm



AC Adapter Jack (Optional)

Allows power supply from an outlet



Last 3 Results Average

Automatically calculates the average of last three readings stored in memory



Latex-free Cuff

Reduces the risk of allergic reaction



WHO Blood Pressure Classification Indicator

Automatically indicates blood pressure classification after the completion of each measurement



120 Memories in 4 Groups

Automatic memory function records the last 120 measurements in 4 groups for efficient and simplified tracking



Low Battery Detection

Visually indicates when the battery has been discharged below a specific level



Auto Shut Off

Shuts down the device when idle for 3 mins, thereby extending the battery life



Deluxe Carry Case

Ensures convenience during travel and storage

Product Specification


Product Description	Arm-type Fully Automatic Digital Blood Pressure Monitor	
Model	BPL 120/180 B1	
Display	LCD Digital Display, Size: 84.1mm X 55.1mm	
Measurement Method	Oscillometric Method	
Measurement Range	Pressure	0mm Hg to 300mm Hg
	Pulse	30 to 180 beats/minute
Measurement Accuracy	Pressure	±8 mmHg (±1.1 kPa) or ±5%
	Pulse	±5%
Pressurization	Automatic Pressurization	
Memory	120 Memories in Four Groups with Date and Time	
Function	Irregular Heartbeat Detection	
	WHO Classification Indicato	
	Last 3 Results Average	
	Low Battery Detection	
	Automatic Power Off	
Power Source	4 AA batteries or Medical AC Adapter (Recommended, not provided)	
Battery Life	Approximately 2 months at 3 tests per day	
Unit Weight	Approx. 405g (Excluding battery)	
Unit Dimensions	Approx. 110 x 155 x 70mm (L x W x H)	
Cuff Circumference	Approx.135 (W)×485(L) mm	
Operating Environment	Temperature	10°C to 40°C
	Humidity	15% to 90% RH
	Pressure	Atmospheric Pressure
Storage Environment	Temperature:	-20°C to 55°C
	Humidity	15% to 90% RH
Classificatio	Internal powered equipment, Type BF, Cuff is the applied par	

Our reach, quick response time
and easy availability of spares
and consumables make us more
accessible anywhere, anytime.

Always close to you

Toll Free: 1800-4252355

From India to the World



CERTIFIED ISO 13485:2003, ISO 9001:2008 COMPANY

BPL Medical Technologies Private Limited

11th KM, Bannerghatta Road,
Arakere, Bangalore - 560076

Ph: +91 80 26484388/ 2648 4350

Toll Free: 1800-4252355

Website: www.bplmedicaltechnologies.com

For Enquiries: sales.medical@bpl.in

Good at Heart
A BPL Medical Technologies Initiative

Wake up to
Wellness
A BPL Medical Technologies Initiative



Happier Living Everyday