

BPL Make Digital Holter 12 Channel & 3 Channel Recorder, Model-TRACK 48 with Holter Software:-

- **LCD Display:-** Provides on demand ECG view/confirmation programming instructions/visual confirmation of patient hook up.
- 12-Lead/3- Channel digital Recording.
- Recording Time: 3CH:upto 96 hrs;12CH up to 48 hrs
- Recording can continue up to 100 days for 3CH and more than one month for 12CH by battery replacement during recording. Recording won't be affected during a battery replacement.
- **Compact Size, Light Weight** 68mm x 53mm x 16mm. 42 grams (excluding Battery)
- Data Transfer: USB or Card Reader
- Comprehensive Holter Software
- Holter software – Easy to use and Understand.

