



Experience the Precision



SURGICAL DIATHERMY

SURGIX-E3 DUAL

Happier Living Everyday



Monopolar Cut

4 Cut Modes: Soft, Pure, Blend & Endo cut



Monopolar Coagulation

3 Coag Modes: Soft, Fulgurate & Spray Coag



Bipolar

2 Bipolar Modes: Precise & Macro Mode



Channels

Monopolar: 2 CH. - 1 Hand & 1 Foot operated
Bipolar: 1 CH. - Foot operated



Smart Bipolar

Audio Feedback after completion of coag to reduce charring & sticking



Advanced Feedback Technology

ESU maintains the set power over wider range of tissue by continuous feedback, no need to adjust the power as tissue changes



Safety Monitoring

Ensures utmost safety with patient plate contact quality monitoring and HF leakage monitoring



Personalized Settings

Allow to personalize settings with upto 10 programs in memory



Self Test

The machine conducts automatic self-test on startup



Application Areas

General Surgery, Gynecology, Laparoscopy, GI Endocut, Papillotomy, Polypectomy, Orthopedic, Plastic surgery, Dermatology, Neurosurgery, ENT, Ophthalmology, Cardiothoracic, Under Water Cutting