

ECUBE8 DIAMOND

Everlasting Productivity

E-CUBE 8 Diamond, the most rational choice for you with excellent imaging performance


Providing both convenience of use and a prompt and accurate diagnostic process, E-CUBE 8 Diamond will become a great partner for your practice. Your patients will receive individual service and care, and the productivity of healthcare will be also improved. E-CUBE 8 Diamond is a smart, reliable product that ensures high performance with outstanding image quality and multi-purpose usability.




Excellent imaging performance

- Flagship-model-grade platform
- Optimal Imaging Suite™ Plus
- Frame rate increasing technology
- Crystal Signature™ convex
- High frequency linear

Frame rate increasing technology



Umbilical Cord in PW



Carotid artery in PW

Simple yet valuable design

- Articulating monitor & Adjustment of CP
- Compact transducer connector
- 21.5-inch full HD LED monitor
- 10.4-inch Touchscreen
- Attached Battery

High resolution transducer



SC1-4HS

SC1-6H


IO8-17

L8-17H


High-performance with clinical versatility




Liver in 2D Mode




Renal vessels in Color Flow Mode



Wrist in 2D Mode



Interventricular septum in Tissue Doppler Imaging



Breast with Elastography



Fetal face in Live HQ™

Transducers

Convex	SC1-4H / SC1-4HS / SC1-6H / C1-6CT C5-8NT
Linear	L3-12H ^{WD} / L3-12T / IO8-17T / L8-17H
Phased	P1-5CT / SP3-8T
Volume	VC1-6T / VE3-10H
Endocavity	EC3-10T / EV3-10T / E3-10
Others	CW 5.0 / CW 2.0

Product Specifications

Weight / Height	55kg / 1470~1560 mm
Width / Depth	532 / 787mm
Monitor	21.5 inch
Touch screen	10.4 inch
Transducer	4ea active
Storage	DICOM 3.0 Connectivity Integrated SSD with HDD 6 USB ports

ALPINION MEDICAL SYSTEMS

www.alpinion.com
Copyright©2019 ALPINION MEDICAL SYSTEMS Co., Ltd. All rights reserved.
Standalone clinical images may have been cropped to better visualize pathology.